

# MIS 100- MEDUSE®

COLUMN IMMOBILIZATION BOARD



The MIS 100 Meduse® Immobilization Board is indicated for use in rescue situations that require the total immobilization of the patient's cervical spine. Meduse® has an anatomical shape, in addition to a radio-transparent structure, which allows its use in X-Ray and Magnetic Resonance exams, as well as a structure that allows for rescues in height and floating in water.

## FEATURES

Polyethylene structure with internal wooden reinforcements  
14 handholds for transport  
12 slots for seat belt adjustments  
Allows coupling with different head immobilizers  
Waterproof, durable and easily asepsis material

## OPTIONAL ACCESSORIES

"Spider" type belt  
Individual belts with velcro closure or quick coupling

## TECHNICAL INFORMATION

Total length: 1830mm  
Width: 440mm  
Thickness: 65mm  
Length of the handgrips: 136mm  
Net Weight: 7.0Kg  
Maximum load capacity: 200Kg

## STANDARD/ CERTIFICACIONES

DIN EN 1865  
ABNT NBR 14.561  
CE Certification  
Anvisa registration nº 80438490002  
AFE Anvisa - 8043849

