


# MPS 110- PANDORA I®

FOLDING STRETCHER



The Pandora I® folding stretcher was developed to facilitate the manual transport of victims in the most diverse rescue situations, through a light, practical and compact equipment.

Pandora I® has a foldable structure that allows it to reduce its dimensions for storage, and can be manufactured in the version with 02 wheels and 02 rubber shoes, or in the version with 04 rubber shoes.

## FEATURES

Structure in hard tubular aluminum  
Base in resistant, washable and waterproof vinyl canvas  
Central folding system for dimensional reduction and storage

## STANDARD ACCESSORIES

02 Quick-connect PVC safety belts

## TECHNICAL INFORMATION

Total open length: 1850mm  
Total folded length: 940mm  
Total width: 460mm  
Total open height: 100mm  
Total folded height: 130mm  
Net Weight: 6.00Kg  
Carrying Capacity: 160.00Kg  
Vinyl Color: Black

## STANDARD/ CERTIFICATIONS

DIN EN 1865 -  
CE Certification  
Anvisa Registration nº 80438490002  
AFE Anvisa - 8043849

